



# **mb-fibreglass**

*Supplying Quality Fibreglass & Moulding Products*

## **PLASTER OF PARIS BANDAGES**

---

### ***Soaking***

Soak one roll bandage together with its spool obliquely into tepid clean water (approx 25C) about 5 seconds, until bubbles stop rising in the water.

---

### ***Squeezing***

Take the bandage gently out of the water and squeeze out the excess water with hands pressing the bandage softly from both edges towards the middle.

---

### ***Wrapping***

Wrap the bandage layer upon layer immediately after squeezing and don't wrap the bandage too tightly.

---

### ***Shaping***

Shape the bandage by hand continuously during wrapping in order to make bandage adhere evenly and remove bubbles. Do not attempt to shape the bandage once it has begun to coagulate.

---

### ***Setting***

The normal setting time is 2-5 minutes, but it can be changed in the following ways :

To accelerate the bandage setting, comparatively warm water may be used, or if necessary a hair dryer may also be used.

The setting time may be retarded approx 4 minutes by adding 2mg sodium sulfite (anhydrous) per 1000ml water before soaking.

---

### ***Storage***

Keep the bandages in well sealed containers, protect from moisture and store in a cool dry place.

*These FAQ's are to be used for general information only. MB Fibreglass will not be held responsible for the accuracy of any information. MB Fibreglass will not be held responsible for any damage caused by the use of these products. Many of the products we supply are industrial grade chemicals and as such it is the customers responsibility to ensure that take all necessary precautions to avoid harm or injury. Always read the technical datasheets and Safety sheets before using any item.*